



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Thermo Scientific TSX Series -20 °C High- Performance Enzyme Freezers without Automatic Defrost

**€9,710.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

TSX Series -20 °C Enzyme Freezers without Auto Defrost Thermo Scientific™ TSX Series -20 °C Enzyme Freezers without Auto Defrost are a sustainable option for enzyme storage. Product Line: TSX Series Certification/Compliance: UL, cUL, CE Capacity (Anglo-American): 29.2 cu. ft. Capacity (metric): 827 L Electrical Requirements: 208-230 V, 50 Hz Frequency: 50 Hz Amperage: 15 A Plug Type: CEE 7/7 Dimensions (L x W x H) Exterior: 38.9 x 34.0 x 78.5 in. (99.0 x 86.4 x 199.4 cm) Dimensions (D x W x H) Interior: 28.5 x 30.0 x 58.0 in. ( 72.3 x 76.2 x 147.3 cm) Weight (Anglo): 484 lb. Shipping Weight (Anglo): 552 lb. Shipping Weight (Metric): 251 kg Display: microprocessor Refrigerant: R290 SmartView Compatibility: Smart-View Compatible Data Recorder: optional Temperature (set point): -25 °C to -15 °C (-20 °C) Interior Depth (Anglo-American): 28.5 in. Exterior height (Anglo-American): 78.5 in. Interior height (Anglo-American): 58.0 in. Inside width (Anglo-American): 30 in. Defrost: Manual Outer width (Anglo-American): 34.0 in. Outside width (metric): 86.4 cm Outside dimensions (D x W x H): 38.9 x 34.0 x 78.5 in. (99.0 x 86.4 x 199.4 cm) Outside depth (Anglo): 38.9 in. Type: High-Performance Freezer Door Type: Fixed Monitoring Options: Wireless No. of Doors: 1 Rack Capacity: 9 Unit Size: Each

## Additional Information

---

No.	LD-5728
Manufacturer (Brand)	Thermo Scientific

