



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Thermo Scientific TSX Series -20 °C- High- Performance Freezers without Automatic Defrost

**€9,840.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

TSX Series -20 °C- high capacity freezers without auto defrost Thermo Scientific™ TSX -20 °C Series high capacity freezers without auto defrost are a sustainable option for clinical and laboratory storage Product Line: TSX Series Certification/Compliance: UL, cUL, CE Capacity (anglo): 29.2 cu. ft. Capacity (metric): 827 L Electrical Requirements: 208-230 V, 50 Hz Frequency: 50 Hz Amperage: 15 A Plug Type: CEE 7/7 Dimensions (D x W x H) Inside: 28.5 x 30.0 x 58.0 in. ( 72.3 x 76.2 x 147.3 cm) Weight (Anglo): 415 lb. Shipping Weight (Anglo): 483 lb. Shipping weight (metric): 219 kg Display: microprocessor Refrigerant: R290 SmartView Compatibility: Smart-Vue Compatible Data Recorder: optional Temperature (set point): -25 °C to -15 °C (-20 °C) Interior Depth (Anglo-American): 28.5 in. Exterior height (Anglo-American): 78.5 in. Interior height (Anglo-American): 58 in. Interior width (Anglo-American): 30 in. Defrost: Manual Outer width (Anglo-American): 34.0 in. Outside width (metric): 86.4 cm Outside dimensions (D x W x H): 38.9 x 34.0 x 78.5 in. (99.0 x 86.4 x 199.4 cm) Outside depth (Anglo): 38.9 in. Type: High-Performance Freezer Eco-friendly Features: Energy efficient Door Type: Fixed Monitoring Options: Wireless No. of Doors: 1 Rack Capacity: 4 Unit Size: Each

## Additional Information

No.	LD-5729
Manufacturer (Brand)	Thermo Scientific

