



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Thermo Scientific TSX Series High-Performance Freezers -30 °C with Automatic Defrost

**€10,678.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

TSX Series -30 °C High-Performance Freezers with Automatic Defrost Thermo Scientific™ TSX Series -30 °C High-Performance Freezers with Automatic Defrost are a sustainable option for clinical and laboratory storage. Product Line: TSX Series Certification/Compliance: UL, cUL, CE Certification/Compliance: cUL Certification/Compliance: UL Capacity (Anglo-American): 51.1 cu. ft. Capacity (metric): 1447 L Electrical Requirements: 208-230 V, 50 Hz Frequency: 50/60 Hz Amperage: 15 A (CB) Plug Type: CEE 7/7 Dimensions (L x W x H) Exterior: 37.9 x 56.5 x 78.5 in. (96.2 x 143.5 x 199.4 cm) Dimensions (D x W x H) Interior: 28.5 x 52.5 x 58 in. (72.4 x 133.4 x 147.3 cm) Weight (anglo): 535 Shipping Weight (anglo): 621 Shipping Weight (metric): 282 Display: Microprocessor Refrigerant: R290 SmartView Compatibility: Smart-View Compatible Data Recorder: Optional Warranty: Warranty varies by country, please contact us for details Temperature (set point): -35 ° to -15 °C (-30 °C) Interior Depth (Anglo-American): 28.5 in. Exterior height (Anglo): 78.5 in. Interior height (Anglo-American): 58 in. Interior width (Anglo-American): 52.5 in. Defrost: Automatic Outer width (anglo): 56.5 in. Outer width (metric): 143.5 in. Outer depth (Anglo-American): 37.8 in. Type: High-Performance Freezer Environmentally friendly features: Energy efficient Door design: 2 solid Monitoring options: Wireless No. of Doors: 2 Rack capacity: 8 Unit Size: Each

## Additional Information

---

No.	LD-5739
Manufacturer (Brand)	Thermo Scientific
VGKL number	107925030

