



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific TSX Series High-Performance Freezers -30 °C with Automatic Defrost

€8,710.00
plus VAT &
Shipping

Product Images

Description

TSX Series High-Performance Freezers -30 °C with Automatic Defrost Thermo Scientific™ TSX Series High-Performance Freezers -30 °C with Automatic Defrost are a sustainable option for clinical and laboratory storage. Product Line: TSX Series Certification/Compliance: UL, cUL, CE Capacity (Anglo): 29.2 cu. ft. Capacity (metric): 827 L Electrical Requirements: 208-230 V, 50 Hz Frequency: 50 Hz Amperage: 15 A Plug Type: CEE 7/7 Dimensions (L x W x H) Exterior: 96.2 x 86.4 x 199.4 cm (37.9 x 34.0 x 78.5 in.) Dimensions (D x W x H) Interior: 28.5 x 30.0 x 58.0 in. (72.3 x 76.2 x 147.3 cm) Weight (Anglo): 371 lb. Shipping Weight (Anglo): 439 lb. Shipping weight (metric): 199 kg Display: microprocessor Refrigerant: R290 SmartView Compatibility: Smart-View Compatible Data Recorder: optional Temperature (set point): -35 ° to -15 °C (-30 °C) Interior Depth (Anglo-American): 28.5 in. Exterior height (Anglo-American): 78.5 in. Inside height (Anglo-American): 58.0 in. Inside width (Anglo-American): 30 in. Defrost: Automatic Outer width (Anglo-American): 34.0 in. Outside width (metric): 86.4 cm Outside dimensions (D x W x H): 37.9 x 34.0 x 78.5 in. (96.2 x 86.4 x 199.4 cm) Outside depth (Anglo): 37.8 in. Type: High-Performance Freezer Eco-friendly Features: Energy efficient Door Type: Fixed Monitoring Options: Wireless No. of Doors: 1 Rack Capacity: 8 Unit Size: Each

Additional Information

No.	LD-5730
Manufacturer (Brand)	Thermo Scientific
VGKL number	107923030

