



Thermo Scientific Variomag Mini Stirrer

Product Images

Description

Variomag Mini Stirrer Thermo Scientific Variomag Mini stirrers are designed for use in the smallest of spaces.
Certification/Compliance: CE Description: Mini w / Standard Tele -Module ; 230 Euro Dimensions (L x W x H): 12 x 12 x 5cm (0.47 x 0.47 x 0.19 in.) Height (anglo): 0.19 in. Length (anglo): 0.47 in. Length (metric): 12 mm Connector type: Euro Protection class: IP68 Temperature range (metric): -10 ° to 56 °C (operating temp) Type: stirrer Warranty and service: 5 years Weight (anglo): 0.04 lb. Width (anglo): 0.47 in. Width (metric): 12 mm Capacity: 5 L Shape: Square Speed: 130 - 1000rpm Stirring Positions: 1 Unit Size: Each

Additional Information

No.	LD-1127
Manufacturer (Brand)	Thermo Scientific
VGKL number	517940163

