



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bühler® roller mixer RM 2

€2,322.00
plus VAT &
Shipping

Product Images



Description

Description Space-saving roller mixer with gentle tilting and rolling motion Scope of delivery Basic unit including 6 removable plastic rollers and collecting tray Technical data Shaking motion: rotating, rocking Shaking surface: 200 x 300 mm Dimension of rollers: 300 x Ø 32 mm Number of rollers: 6 pcs Max. Load: 2 kg Speed: 5 - 80 min⁻¹ (adjustable in steps of 5) Rocking angle: 3° fixed Stroke at roller end: 15 mm Running time: 0 - 99 h 59 min LED display: for speed and running time Permissible humidity: ~ 85% non-condensing Ambient temperature: 5°C - 45°C Electrical connection: 100 - 240 V, 50/60 Hz Protection class: IP 21 Heat output: approx. 5 W External dimensions (W x D x H): 244 x 465 x 115 mm Weight: 7.3 kg Features: Suitable for all common tube sizes, bottles and other universal vessels up to 80 mm in diameter due to removable rollers

Additional Information

No.	NB-0134
Manufacturer (Brand)	Bühler
VGKL number	106168000
Transport temperature	Room temperature

