



Krüss® Binocular Microscope MBL2000-LED

Product Images

Description

With 7 LED lamp and brightness control 3.6W power LED similar to 50W halogen bulb with 45° inclined view Eyepieces: 10x plano eyepieces Objectives: achromatic 4x/NA 0.10; 10x/NA 0.25; 40x /0.65; 100x/1.25 oil immersion. Size stage: 130x140 mm Movement stage: 77x50 mm Illumination: 6 V, 20 W Power supply: 240 V, 20 W, 50/60 Hz Double lens Abbe condenser N.A. 1.25 with iris diaphragm, swing-in filter holder, height adjustment. Glass filters: blue, green, yellow

Additional Information

No.	KS-0350
Manufacturer (Brand)	Krüss
VGKL number	106110108

