



Krüss® Binocular Microscope MBL2000-PL

Product Images



Description

with 45° inclined view Eyepieces: 10x plano eyepieces Objectives: planachromatic 4x/0.10; 10x/0.25; 40x/0.65; 100x/1.25 oil immersion. Sturdy metal housing with coaxial fine and coarse drive with adjustable end stop. Stage size: 130x140 mm Stage movement: 77x50 mm Built-in illumination 6 V 20 W Double lens Abbe condenser NA 1.25 with iris diaphragm and swing-in filter holder. Height adjustment. Glass filters: blue, green, yellow Power supply 230 V, 50/60 Hz

Additional Information

No.	KS-0351
Manufacturer (Brand)	Krüss

