



Krüss® Trinocular Microscope MBL2000-T-LED

Product Images

Description

With 45° inclined view and additional photo tube. Eyepieces: 10x plano eyepiece pair Objectives: achromatic 4x/0.10; 10x/0.25; 40x/0.65; 100x/1.25 oil immersion. Sturdy metal housing with coaxial fine and coarse drive with adjustable end stop. Size object stage: 130x140 mm Movement object stage: 77x50 mm Built-in LED illumination brightness control. Double lens Abbe condenser N.A. 1.25, with iris diaphragm, swing-in filter holder. Filters: blue, green, yellow. Power supply 230 V 50/60 Hz

Additional Information

No.	KS-0357
Manufacturer (Brand)	Krüss
VGKL number	106110110

