



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



WELCH® Suction/Filtration System 2511, 9.2, l/min, 292 mbar, 230 V, 50/60 Hz

€815.00
plus VAT &
Shipping

Product Images



Description

System includes: autoclavable pressure vessel, adapters and tubing, in-line filter hydrophobic, potentiometer, pressure gauge and automatic shut-off valve. Pumping speed 50/60 Hz 9.2/11 l/min Suction rate for liquids 15 ml/s Noise level 45 dB(A) Collection container plastic 1.2 l Option foot switch, handpiece with pipetting set

Additional Information

No.	WL-0808
Manufacturer (Brand)	WELCH
VGKL number	109242511
Transport temperature	Room temperature

