



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle® Dosing pump Dosi-Pump, 100 ml/stroke, fixed

€86.60 As low

as

€65.60

**plus VAT &
Shipping**



Description

Standard dosing approx. 100 ml per stroke, the pump has a fixed outlet. Dosing pumps suitable for drums, canisters, balloons or bottles, for dosing liquids of different viscosity. The dosing pumps are mounted by you in your existing container closure (see table mounting hole). For this purpose, an appropriate hole is drilled in the closure, the pump is pushed through and screwed against it from below. A suitable container connection with location hole is not included in the scope of delivery. Wetted parts made of PP, PE and stainless steel

Additional Information

No.	KL-1101
Manufacturer (Brand)	Bürkle
VGKL number	100457101
Transport temperature	Room temperature

