




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bürkle® Dosing pump Dosi-Pump, 30 ml/stroke, fixed

**€27.40** As low

as

**€20.75**

**plus VAT &  
Shipping**



## Description

---

Standard dosing approx. 30 ml per stroke, the pump has a fixed outlet. Dosing pumps suitable for drums, canisters, balloons or bottles, for dosing liquids of different viscosity. The dosing pumps are mounted by you in your existing container closure (see table mounting hole). For this purpose, an appropriate hole is drilled in the closure, the pump is inserted through it and screwed against it from below. A suitable container connection with location hole is not included in the scope of delivery. Wetted parts made of PP, PE and stainless steel

## Additional Information

---

No.	KL-1100
Manufacturer (Brand)	Bürkle
VGKL number	100457031
Transport temperature	Room temperature

