



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle® Dosing pump Dosi-Pump, 4ml/stroke, immersion depth 25cm

€10.10 As low as
€7.55
**plus VAT &
Shipping**

Product Images



Description

Standard dosing approx. 4 ml per stroke. Dosing pumps suitable for drums, canisters, balloons or bottles, for dosing liquids of different viscosity. The dosing pumps are mounted by you in your existing container closure (see table mounting hole). For this purpose, an appropriate hole is drilled in the closure, the pump is inserted through it and screwed against it from below. A suitable container connection with location hole is not included in the scope of delivery. Wetted parts made of PP, PE and stainless steel

Additional Information

No.	KL-1098
Manufacturer (Brand)	Bürkle
VGKL number	100457004
Transport temperature	Room temperature

