




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




---

## KNF® Dosing Pump SIMDOS 02, FEM1.02FT.18S PTFE

**€2,938.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

Chemically resistant diaphragm metering pump SIMDOS 02 FT.18 S - Delivery rate from 0.03 - 20 ml/min / Pressure head 60 mWS, suction head 2 mWS - Pump head in material version PTFE - Diaphragm FFKM, valves made of FFKM - Delivery rate manually adjustable (version S) - Additional safety diaphragm for maximum protection - Easy replacement of the delivery diaphragm by activating the maintenance command stored in the operating program

Technical data Application Dosing of liquids Delivery rate (ml/min) with water of 20 °C at demand head zero 0,03 - 20 Pressure head (mWS) 60 Suction head (mWS) 2 Hose connection (mm) ID1.6 / AD3.2 Permissible media and ambient temperature Ambient temperature +5 ... +40 °C Media temperature +5 ... +80 °C Weight (kg) 0.9 Dimensions W x H x D (mm) 93 x 144 x150 Material pump head PTFE Diaphragm FFKM Valves FFKM Included in delivery: -Dosing pump - Suction filter -2 meter connection hose -Power supply Accessories (not included in delivery) Stand bracket Order no. 160474 Wall bracket Order no. 160473 In-line filter FS 60 T PVDF Mesh size 70 µm Order no. 165210 FS 60 X PEEK Mesh size 35 µm Order no. 165212 Pressure retaining valve JFDV30FTZ: 166288

## Additional Information

No.	KF-0250
Manufacturer (Brand)	KNF
VGKL number	106271283
Transport temperature	Room temperature

