




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KNF® Dosing Pump SIMDOS 10 FEM, .10KT.18S PP

**€1,600.00**  
**plus VAT &**  
**Shipping**



## Description

---

Chemically resistant diaphragm metering pump SIMDOS 10 KT.18 S2 - Delivery rate from 1 - 100 ml/min / Pressure head 60 mWS, suction head 3 mWS - Pump head in material version PP - Diaphragm made of PTFE, Valves made of FFKM - Delivery rate manually adjustable (version S) - Additional safety diaphragm for maximum protection - Easy replacement of the delivery diaphragm by activating the maintenance command stored in the operating program

Technical data Application Dosing of liquids Delivery rate (ml/min) with water of 20 °C at zero delivery head 1 - 100 Pressure head (mWS) 60 Suction head (mWS) 3 Hose connection (mm) ID4 / AD6 Permissible media and ambient temperature Ambient temperature +5 ... +40 °C Media temperature +5 ... +80 °C Weight (kg) 0.9 Dimensions W x H x D (mm) 93 x 144 x150 Material Pump head PP Diaphragm PTFE Valves FFKM Included in delivery: -Dosing pump -Suction filter -2 meter connection hose -Power supply Accessories (not included in scope of delivery) Stand bracket Order no. 160474 Wall bracket Order no. 160473 In-line filter FS 25 T PVDF Mesh size 70 µm Order no. 165211 FS 25 X PEEK Mesh size 35 µm Order no. 165213

## Additional Information

---

No.	KF-0259
Manufacturer (Brand)	KNF
VGKL number	106271181
Transport temperature	Room temperature

