



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KNF® Dosing Pump SIMDOS 10, FEM1.10TT.18RCP PVDF

€2,233.00
plus VAT &
Shipping



Description

Chemically resistant diaphragm metering pump SIMDOS 10 TT.18 RC P2 - delivery rate from 1 - 100 ml/min / pressure head 60 mWS, suction head 3 mWS - pump head in material version PVDF - diaphragm made of PTFE, Valves made of FFKM - Delivery rate adjustable manually and by means of external control as well as including interface RS 232 (version RC Plus) - Additional safety diaphragm for maximum protection - Easy replacement of the delivery diaphragm by activating the maintenance command stored in the operating program Technical data Application Dosing of liquids Delivery rate (ml/min) with water of 20 °C at delivery head zero 1 - 100 Pressure head (mWS) 60 Suction head (mWS) 3 Hose connection (mm) ID4 / AD6 Permissible media and ambient temperature Ambient temperature +5 ... +40 °C Media temperature +5 ... +80 °C Weight (kg) 0.9 Dimensions W x H x D (mm) 93 x 144 x150 Material Pump head PVDF Diaphragm PTFE Valves FFKM Included in delivery: -Dosing pump -Suction filter -2 meter connection hose -Power supply Accessories (not included in scope of delivery) Stand bracket Order no. 160474 Wall bracket Order no. 160473 Foot slider for RC version (RC = manual operation and external control) Order no. 155872 In-line filter FS 25 T PVDF Mesh size 70 µm Order no. 165211 FS 25 X PEEK Mesh size 35 µm Order no. 165213 Hose (material PTFE): 019241

Additional Information

No.	KF-0258
Manufacturer (Brand)	KNF
VGKL number	106271185
Transport temperature	Room temperature

