



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bohlender® BOLA cordless acid and alkali pumps approx. 6 ltr/min, L 400 mm

€418.03
plus VAT &
Shipping

Product Images



Description

BOLA Cordless Acid and Lye Pumps Made of polypropylene, PTFE, Hastelloy®, driven by two standard 1.5 volt batteries (rechargeable batteries are recommended). Dimensions: Suction tube length mm 400 Suction tube Ø, mm approx. 25

Additional Information

No.	ER-1550
Manufacturer (Brand)	Bohlender
VGKL number	291787001
Transport temperature	Room temperature

