



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle® Electr. barrel pump Petro, with automatic nozzle

€674.05
plus VAT &
Shipping

Product Images



Description

Electric barrel pump for pumping diesel, EL heating oil, undiluted radiator antifreeze and other liquids of hazard class AIII (VbF). The integrated siphon protection prevents leakage in case of a defective hose. Available for 230 V. Optionally with standard or self-switching automatic nozzle. Electronic flow meter with totalizer and calibration option (uncalibrated) available as accessory. With manual suction pump VDE-tested, motor with thermal switch fuse Barrel thread M 64x4 and 2" suction hose with filter 1.60 m long Dispensing hose made of Perbunan, 4 m long, with anti-kink spring and nozzle Not suitable for watery/diluted liquids!

Additional Information

No.	KL-1144
Manufacturer (Brand)	Bürkle
Transport temperature	Room temperature

