



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle® Hand lever pump Alu, electr. conductive, M64x4 / R2"

€228.25
plus VAT &
Shipping

Product Images



Description

Lever pump for petrochemicals made of sturdy die-cast aluminum, valves and piston made of plastic (PA). Suitable for barrel thread M 64x4 or R 2" With locking device (for padlock) Immersion tube length 84 cm Delivery includes hose set, hose length 1.5 m, electrically conductive, with discharge elbow and anti-kink spring

Additional Information

No.	KL-1123
Manufacturer (Brand)	Bürkle
Transport temperature	Room temperature

