



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle® Lever Pump Compact, 500ml/stroke, 2" BSP, 100 cm

€83.90
plus VAT &
Shipping

Product Images



Description

Compact, self-priming lever pump made of polypropylene, suitable for light acids and alkalis. The telescopic tube of the lever pump is continuously extendable up to 100 cm, it can thus be used for 50 - 200 l drums. Polypropylene, seals made of NBR Delivery rate up to approx. 500 ml per stroke Outlet bend can be rotated through 360° With integrated barrel screw connection for G2" thread Attention! "List of chemical resistances" must be observed!

Additional Information

No.	KL-1107
Manufacturer (Brand)	Bürkle
Transport temperature	Room temperature

