



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bürkle® lever pump, blue, galvanized piston rod

**€40.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The lever pump is made of polypropylene with a galvanized steel piston rod and FKM seals. This manual lever pump can be used to remove quantities as needed from drums without much effort. Ideal for occasional use. Applications of the lever pump: water, lubricants, light oils, as well as water-based chemicals and light acids and alkaline solutions. With integrated barrel screw connection for 1 1/4" and 2" fine thread Delivery rate 350 ml per stroke Immersion depth infinitely variable from 65-90 cm "Note list of chemical resistances"!

## Additional Information

---

No.	KL-1109
Manufacturer (Brand)	Bürkle
Transport temperature	Room temperature

