



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle® Laboratory pump with mains operation, PP/PVC, 70 cm, 75l/min

€1,056.35
plus VAT &
Shipping

Product Images



Description

Lightweight and powerful pump for filling aqueous media up to 300 mPas from canisters, drums and narrow-necked containers from 32 mm opening diameter. PP Splash-proof motor 230 V/200 W (IP44) Integrated safety button as on/off switch Sealless pumping mechanism for acids and alkalis Threaded connection 1" with hose nozzle 3/4" Incl. 1.5 m PVC hose 3/4", PP filler and suspension device Attention! "List of chemical resistances" must be observed!

Additional Information

No.	KL-1128
Manufacturer (Brand)	Bürkle
Transport temperature	Room temperature

