



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bürkle® PumpMaster aqueous liquids, PP/EPDM, red

**€178.65** As low  
as  
**€150.65**  
plus VAT &  
Shipping

## Product Images

---



## Description

---

PumpMaster, canister and barrel pump for non-aggressive, aqueous liquids. The manual pressure pump is ideal for filling, transferring and emptying canisters and drums. An overpressure is created by means of a hand piston, which forces the liquid out of the container. The discharge tap also allows small quantities of liquid to be removed. The PumpMaster canister and barrel pump is mounted on the canister with a flexible plug that can be variably adjusted in diameter. No special thread adapters are required. PumpMaster fits canister openings from 46.5 mm to 60 mm in diameter. The immersion tube has a length of 95 cm and can be shortened depending on the canister size. Complete with three plugs in different sizes Not suitable for hydrochloric acid! "Observe list of chemical resistances"!

## Additional Information

---

No.	KL-0573
Manufacturer (Brand)	Bürkle
VGKL number	100452200
Transport temperature	Room temperature

