




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



WELCH® vacuum/pressure pump 2581, 83 l/min, 6.7, mbar

€2,085.00
plus VAT &
Shipping

Product Images



Description

The 2581 piston pump is ideal for pumping steam, and also for applications with vacuum requirements up to 6.7 mbar. The pump provides a continuous, reliable and high vacuum. It is lightweight and compact. The pump's oil-free operation requires little maintenance and comes with a liquid trap and mounted vacuum gauge. Note: Not recommended for organic, acidic and basic vapors.

Additional Information

No.	WL-0834
Manufacturer (Brand)	WELCH
VGKL number	109242581
Transport temperature	Room temperature

