



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KNF® Diaphragm Vacuum Pump N842.3FT.18

€3,396.00
plus VAT &
Shipping

Product Images



Description

Diaphragm vacuum pump N 842.3 FT.18 - Delivery rate 2.04 m³/h / Ultimate vacuum 2 mbar abs. - High steam and condensate compatibility - PTFE pump head in combination with PTFE-coated diaphragm are ideal for very for aggressive/corrosive gases and vapors Technical data Application Rotary evaporation, distillation Delivery rate (m³/h) at atm. Pressure 2.04 Ultimate vacuum (mbar abs.) 2 Operating overpressure (bar) 1 Hose connections (mm) ID10 Permissible media and ambient temperature +5 ... +40 °C Weight (kg) 13.4 Dimensions W x H x D (mm) 167 x 223 x 341 Material Pump head PTFE Diaphragm PTFE-coated Valves FFPM Accessories (optional): Vacuum chemistry hose 028187 Noise damper 056366 (mounting by means of hose on outlet side)

Additional Information

No.	KF-0158
Manufacturer (Brand)	KNF
VGKL number	566240260
Transport temperature	Room temperature

