




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KNF® Diaphragm Vacuum Pump N860.3FT.40.18

€4,320.00
plus VAT &
Shipping

Product Images



Description

Diaphragm vacuum pump N 860.3 FT.40.18 - Delivery rate 3.6 m³/h / Ultimate vacuum 4 mbar abs. - Integrated KNF self-drying system ensures rapid removal of condensate from the pump heads without changing the vacuum. This significantly reduces process time and protects the pump heads. - Chemically resistant and therefore also suitable for use with very aggressive/corrosive gases and vapors Technical data Application Rotary evaporation, distillation, vacuum furnace, vacuum concentrator Delivery rate (m³/h) at atm. pressure 3.6 Pressure 3.6 Ultimate vacuum (mbar abs.) 4 Operating overpressure (bar) 1 Hose connections (mm) ID12 Permissible media and ambient temperature +5 ... +40 °C Weight (kg) 14.8 Dimensions W x H x D (mm) 291 x 278 x 331 Material Pump head PTFE Diaphragm PTFE-coated Valves FFPM

Additional Information

No.	KF-0147
Manufacturer (Brand)	KNF
VGKL number	106278603
Transport temperature	Room temperature

