



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KNF® Vacuum pump LABOPORT N938.50KT.18

€1,453.00
plus VAT &
Shipping

Product Images



Description

Mini diaphragm vacuum pump N 938.50 KT.18 - Delivery rate 1.8 m³/h / Ultimate vacuum 15 mbar abs. - Parallel and in-line connection of the two pump heads ensures very fast evacuation - PTFE-coated diaphragm for aggressive/corrosive gases and vapors Technical data Application Filtration, SPE, degassing Delivery rate (m³/h) at atm. pressure 1.8 Ultimate vacuum (mbar abs.) 15 Pressure 1.8 Ultimate vacuum (mbar abs.) 15 Operating overpressure (bar) 0.5 Hose connections (mm) ID10 Permissible media and ambient temperature +5 ... +40 °C Weight (kg) 6.8 Dimensions W x H x D (mm) 110 x 212 x 317 Material pump head PPS Diaphragm PTFE-coated Valves FFPM Accessories Noise damper Order no. 007006 Hose nipple G1/8 ID10 PVDF Order no. 112004 Fine adjustment valve with vacuum gauge, suction side for adjusting vacuum Order no. 112432 Accessories (optional): Vacuum chemistry hose 028187 Silencer 056366 (mounting by means of hose on outlet side)

Additional Information

No.	KF-0246
Manufacturer (Brand)	KNF
VGKL number	106270150
Transport temperature	Room temperature

