



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KNF® Vacuum System SC 820

€5,041.00
plus VAT &
Shipping

Product Images



Description

Vacuum pumping system SC 820 - Delivery rate 1.2 m³/h / Ultimate vacuum 8 mbar abs. - Vacuum system consisting of chemically resistant diaphragm vacuum pump, base plate, condenser, separator and vacuum controller Technical data Application Rotary evaporation Delivery rate (m³/h) at atm. Pressure 1.2 Ultimate vacuum (mbar abs.) 8 Operating overpressure (bar) 1 Hose connections (mm) pneumatic: ID10 Coolant: ID8 Permissible media and ambient temperature +5 ... +40 °C Weight (kg) 16.0 Dimensions W x H x D (mm) 289 x 506 x 397 Material Pump head PTFE Diaphragm PTFE-coated Valves FFPM Accessories Coolant valve - G1/2, ID8 Order no. 045075 Accessories (optional) Vacuum chemistry hose ID 10mm (delivered in whole meters): 028187 Vacuum chemistry hose ID 6mm (delivered in whole meters): 055535 Small flange (in connection with vacuum hose): 045419

Additional Information

No.	KF-0138
Manufacturer (Brand)	KNF
VGKL number	566240620
Transport temperature	Room temperature

