



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



VACUUBRAND® Chemistry Diaphragm Pump MD 12C NT +EK, three-stage, 230 V / 50-60 Hz, CEE power cable

€8,730.00
plus VAT &
Shipping

Product Images



Description

Chemistry vacuum system MD 12C NT +EK --- --- Max. Pumping speed 50/60 Hz (m³/h): 11.1/12.3 --- Ultimate vacuum (abs.) (mbar): 2 --- Ultimate vacuum (abs.) with gas ballast (mbar): 4 --- Connection suction side: small flange KF DN 25 --- Connection discharge side: hose shaft DN 10 mm --- Connection cooling water: 2 x hose shaft DN 6-8 mm --- Dimensions (L x W x H), approx. (mm): 528x387x395 --- Weight, approx. (kg): 29.1 --- Rated voltage / mains frequency: 230 V/50-60 Hz --- Plug: CEE --- Motor rated power (kW): 0.44 --- Scope of delivery: Pumping unit completely assembled, ready for operation, with instructions. Old item number: 743500

Additional Information

| | |
|-----------------------|------------------|
| No. | VZ-1733 |
| Manufacturer (Brand) | VACUUBRAND |
| Transport temperature | Room temperature |

