



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



VACUUBRAND® chemistry diaphragm pump MD 4C EX three-stage, 230 V / 50 Hz cable without plug

€8,730.00
plus VAT &
Shipping

Product Images



Description

Chemistry diaphragm pump MD 4C EX --- --- Max. Pumping speed 50/60 Hz (m³/h): 3.7/ --- Ultimate vacuum (abs.) (mbar): 3 --- Ultimate vacuum (abs.) with gas ballast (mbar): 10 --- Connection suction side: Small flange KF DN 25 --- Connection discharge side: Small flange KF DN 16 --- Dimensions (L x W x H), approx. (mm): 440x265x305 --- Weight, approx. (kg): 29.3 --- Nominal voltage / mains frequency: 230 V/50 Hz --- Connector: Open strands --- Motor nominal power (kW): 0.25 --- Scope of delivery: Pump completely assembled, ready for operation, with instructions, with 2 m cable with open strands Old item number: 696930

Additional Information

No.	VZ-1334
Manufacturer (Brand)	VACUUBRAND
Transport temperature	Room temperature

