



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



VACUUBRAND® Chemistry Diaphragm Pump MD 4C NT, three-stage Certification (NRTL): C/US 230 V / 50-

€4,165.00
plus VAT &
Shipping

Product Images



Description

Chemistry diaphragm pump MD 4C NT --- --- Max. Pumping speed 50/60 Hz (m³/h): 3.4/3.8 --- Ultimate vacuum (abs.) (mbar): 1.5 --- Ultimate vacuum (abs.) with gas ballast (mbar): 3 --- Connection suction side: hose shaft DN 10 mm --- Connection discharge side: hose shaft DN 10 mm --- Dimensions (L x W x H), approx. (mm): 325x243x198 --- Weight, approx. (kg): 14.3 --- Rated voltage / mains frequency: 230 V/50-60 Hz --- Plug: CEE --- Motor rated power (kW): 0.25 --- Scope of delivery: pump completely assembled, ready for operation, with instructions. Old item number: 736400

Additional Information

No.	VZ-1611
Manufacturer (Brand)	VACUUBRAND
Transport temperature	Room temperature

