



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



VACUUBRAND® Chemistry Diaphragm Pump ME 16C NT +EK, single stage, 230 V / 50-60 Hz, CEE power cable

€8,650.00
plus VAT &
Shipping

Product Images



Description

ME 16C NT +EK chemistry vacuum system --- --- Max. Pumping speed 50/60 Hz (m³/h): 16.3/18.4 --- Ultimate vacuum (abs.) (mbar): 70 --- Ultimate vacuum (abs.) with gas ballast (mbar): 100 --- Connection suction side: small flange KF DN 25 --- Connection discharge side: hose shaft DN 10 mm --- Connection cooling water: 2 x hose shaft DN 6-8 mm --- Dimensions (L x W x H), approx. (mm): 528x387x395 --- Weight, approx. (kg): 29.1 --- Rated voltage / mains frequency: 230 V/50-60 Hz --- Plug: CEE --- Motor rated power (kW): 0.44 --- Scope of delivery: pump completely assembled, ready for operation, with instructions. Old item number: 741500

Additional Information

No.	VZ-1711
Manufacturer (Brand)	VACUUBRAND
Transport temperature	Room temperature

