



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



VACUUBRAND® Diaphragm pump ME 4R NT, single-stage, 230 V / 50-60 Hz, CEE power cable Certification

€1,925.00
plus VAT &
Shipping

Product Images



Description

Diaphragm pump ME 4R NT --- --- Max. Pumping speed 50/60 Hz (m³/h): 3.8/4.2 --- Ultimate vacuum (abs.) (mbar): 100 --- Vacuum controller: pressure gauge --- Connection suction side: hose shaft DN 10 mm --- Connection discharge side: hose shaft DN 10 mm --- Dimensions (L x W x H), approx. (mm): 243x239x290 --- Weight, approx. (kg): 11.5 --- Nominal voltage / mains frequency: 230 V/50-60 Hz --- Plug: CEE --- Motor nominal power (kW): 0.18 --- Scope of delivery: pump completely assembled, ready for operation, with instructions. Old item number: 731100

Additional Information

No.	VZ-1490
Manufacturer (Brand)	VACUUBRAND
Transport temperature	Room temperature

