



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**VACUUBRAND® diaphragm pump MZ 2D NT,  
two-stage, certification (NRTL): C/US 230 V /  
50-60 Hz,**

**€1,965.00  
plus VAT &  
Shipping**

## Product Images



## Description

---

Diaphragm pump MZ 2D NT --- --- Max. Pumping speed 50/60 Hz (m<sup>3</sup>/h): 2.3/2.5 --- Ultimate vacuum (abs.) (mbar): 4 --- Connection suction side: small flange KF DN 16 --- Connection discharge side: silencer / G1/4" --- Dimensions (L x W x H), approx. (mm): 243x242x198 --- Weight, approx. (kg): 11.4 --- Nominal voltage / mains frequency: 230 V/50-60 Hz --- Plug: CEE --- Motor nominal power (kW): 0.18 --- Scope of delivery: pump completely assembled, ready for operation, with instructions. Old item number: 732200

## Additional Information

---

No.	VZ-1505
Manufacturer (Brand)	VACUUBRAND
Transport temperature	Room temperature

