



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## VACUUBRAND® VARIO Chemistry Pumping Unit PC 3001 VARIO select EK Peltronic

**€7,620.00**  
**plus VAT &**  
**Shipping**



## Description

---

--- VACUU-SELECT vacuum controller with graphical user interface and predefined vacuum processes simplifies work in the laboratory --- automatic boiling point detection and adjustment of vacuum for short process times --- extremely compact and powerful, very good ultimate vacuum even with open gas ballast valve --- very quiet and low-vibration operation --- excellent environmental compatibility due to low energy consumption and efficient solvent recovery --- Max. Pumping speed 50/60 Hz (m<sup>3</sup>/h): 2.0/ --- Ultimate vacuum (abs.) (mbar): 2.0 --- Ultimate vacuum (abs.) with gas ballast (mbar): 4 --- Vacuum controller: VACUU-SELECT --- Connection suction side: hose shaft DN 6-10 mm --- Connection discharge side: hose shaft DN 8-10 mm --- Dimensions (L x W x H), approx. (mm): 300x370x400 --- Weight, approx. (kg): 11.8 --- Rated voltage / mains frequency: 200-230 V/50-60 Hz --- Plug: IEC plug EN 60320 --- Motor rated power (kW): 0.16 --- Scope of delivery: Pump completely assembled, ready for operation, with instructions.

## Additional Information

---

No.	VZ-1442
Manufacturer (Brand)	VACUUBRAND
Transport temperature	Room temperature

