



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



VACUUBRAND® VARIO Chemistry Pumping Unit PC 3001 VARIO select TE

€5,820.00
plus VAT &
Shipping



Description

--- Vacuum controller VACUU-SELECT with graphical user interface and predefined vacuum processes simplifies work in the laboratory --- automatic boiling point detection and adjustment of vacuum for short process times --- extremely compact and powerful, very good ultimate vacuum even with open gas ballast valve --- very quiet and low-vibration operation --- excellent environmental compatibility due to low energy consumption and efficient solvent recovery --- Max. Pumping speed 50/60 Hz (m³ /h): 2.0/ --- Ultimate vacuum (abs.) (mbar): 2.0 --- Ultimate vacuum (abs.) with gas ballast (mbar): 4 --- Vacuum controller: VACUU-SELECT --- Connection suction side: hose shaft DN 6-10 mm --- Connection discharge side: hose shaft DN 8-10 mm --- Dimensions (L x W x H), approx. (mm): 301x341x493 --- Weight, approx. (kg): 8.7 --- Rated voltage / mains frequency: 200-230 V/50-60 Hz --- Plug: CEE --- Motor rated power (kW): 0.16 --- Scope of delivery: Pump completely assembled, ready for operation, with instructions.

Additional Information

No.	VZ-1440
Manufacturer (Brand)	VACUUBRAND
Transport temperature	Room temperature

