



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



WELCH® Diaphragm pump MP 1201 E 230V 50/60Hz

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Pumping speed 50/60 Hz 138/151 l/min Motor power 370 W Suction port small flange DN 16 KF Pressure port hose shaft DN 8 for hose inner diameter 8 mm Option gas ballast valve, hose shaft DN 8

Additional Information

No.	WL-0222
Manufacturer (Brand)	WELCH
VGKL number	569240392
Transport temperature	Room temperature

