



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



WELCH® Diaphragm pump chemically resistant MPC 901 Z 230V, 50/60Hz

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Pumping speed 50/60 Hz 113/125 l/min Motor power 370 W Suction connection small flange DN 16 KF (hose shaft DN 8 enclosed) Pressure connection hose shaft DN 8 for internal hose diameter 8 mm

Additional Information

No.	WL-0270
Manufacturer (Brand)	WELCH
VGKL number	569240522
Transport temperature	Room temperature

