



Lovibond® VIS Spectrophotometer XD7000

Product Images



Description

The instruments of the XD series are ideally suited for routine and spectral analyses. As true all-rounders, the XD7000 and XD7500 simplify the work routine, especially with automatic method recognition using barcodes, automatic cuvette recognition and more than 150 pre-programmed methods, and also support Good Laboratory Practice (GLP) in addition to analytical quality control. The instruments are available together with the barcoded cuvettes and extensive accessories from a single source - which is also easy on your wallet. High-quality reference beam optics, user-friendliness, easy handling and multilingual global usability round off the versatility of the VIS and UV/VIS spectrophotometers, which are also equipped for flexible on-site use.

Additional Information

No.	CV-0472
Manufacturer (Brand)	Lovibond
VGKL number	107717000

