



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



WLD-TEC® SteriMax basic 230V, infrared inoculation loop sterilizer

€737.00
plus VAT &
Shipping

Product Images



Description

SteriMax basic, 230V - Infrared inoculation loop sterilizer High-speed annealing and sterilization without gas flame. Best suited for all laboratories, anaerobic work environments and safety cabinets! Specially focused infrared light creates an IR HotSpot and sterilizes your inoculating loop at a temperature of 750 to 1000°C in just 5 to 10 seconds with power from the wall socket, instantly and without any preheating! - Immediate readiness for use without heating phase - No risk of burns on the housing (Cool-Touch) - Convenient touch operation - Efficient use of energy through Thermocontrol - Annealing tube made of special quartz glass - Easy cleaning of device and annealing tube

Additional Information

No.	WT-0024
Manufacturer (Brand)	WLD-TEC
VGKL number	103074110
Transport temperature	Room temperature

