



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bandelin® SONOREX DIGIPLUS DL 514 BH Ultrasonic bath with heater 18.7 liters

€2,430.00
plus VAT &
Shipping

Product Images



Description

High-performance ultrasonic bath with power adjustment from 20 to 100 % in 10 % steps and heating for use in the laboratory and process engineering. Possible applications: - Cleaning of sensitive materials - Mixing of samples - Homogenization of samples for analyses - Rapid degassing of solvents for HPLC The proven and particularly robust SONOREX ultrasonic baths clean thoroughly and gently by means of intensive cavitation bubbles even in the hard-to-reach places of the cleaning objects. A user-friendly and easy-to-wipe button control is used for simple setting of the cleaning time and bath temperature. The remaining cleaning time as well as the current set and actual values of the bath temperature can be clearly read on the display. The fully automatic SweepTec function prevents standing waves in the ultrasonic field and achieves particularly uniform cleaning. The degassing function (Degas), which can be activated at the beginning of the cleaning process, additionally enhances the cleaning effect of the ultrasound. Suitable baskets, lids and cleaning and disinfection concentrates are available as accessories for optimally coordinated ultrasonic cleaning. Internal dimensions L x W x D: 325 x 300 x 200 mm Mains connection: 230 V EU plug

Additional Information

No.	BL-0344
Manufacturer (Brand)	Bandelin
VGKL number	100737185
Transport temperature	Room temperature

