



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bandelin® SONOREX DIGITEC DT 1028 F Ultrasonic bath 9.5 liters

**Price on
request
plus VAT &
Shipping**

Product Images



Description

High-performance ultrasonic bath with flat oscillating tank for rapid degassing and sample preparation in the laboratory. Uniform sonication of all samples regardless of their size and arrangement. Homogenization or rapid degassing of samples. User-friendly and easy to wipe button operation is used for easy setting of sonication duration. Holding clamps prevent the laboratory flasks from tilting or floating. The remaining sonication time can be clearly read on the display. The fully automatic SweepTec function prevents standing waves in the ultrasonic field and ensures particularly uniform sonication. Internal dimensions L x W x D: 500 x 300 x 65 mm Mains connection: 230 V EU plug

Additional Information

No.	BL-0159
Manufacturer (Brand)	Bandelin
VGKL number	100733243
Transport temperature	Room temperature

