



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bandelin® SONOREX DIGITEC DT 1028 Ultrasonic bath 28 liters

€2,662.00
plus VAT &
Shipping

Product Images



Description

High-performance ultrasonic cleaners for use in the practice, laboratory and workshop. The proven and particularly robust SONOREX ultrasonic baths clean thoroughly and gently by means of intensive cavitation bubbles even in the hard-to-reach places of the cleaning objects. A user-friendly and easy-to-wipe button operation is used for easy setting of the cleaning time. The remaining cleaning time can be clearly read on the display. The fully automatic SweepTec function prevents standing waves in the ultrasonic field and achieves particularly even cleaning. The degassing function (Degas), which can be activated at the beginning of the cleaning process, further enhances the cleaning effect of the ultrasound. Suitable baskets, lids and cleaning and disinfection concentrates are available as accessories for optimally coordinated ultrasonic cleaning. Internal dimensions L x W x D: 500 x 300 x 200 mm Mains connection: 230 V EU plug

Additional Information

No.	BL-0166
Manufacturer (Brand)	Bandelin
Transport temperature	Room temperature

