



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bandelin® SONOREX DIGITEC DT 510 F Ultrasonic bath 4.3 liters

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

Flat ultrasonic bath for rapid degassing and sample preparation in the laboratory. Uniform sonication of all samples regardless of their size and arrangement. Homogenization or rapid degassing of samples. User-friendly and easy to wipe button operation is used for easy setting of sonication duration. Holding clamps prevent the laboratory flasks from tilting or floating. The remaining sonication time can be clearly read on the display. The fully automatic SweepTec function prevents standing waves in the ultrasonic field and ensures particularly uniform sonication. Internal dimensions L x W x D: 300 x 240 x 65 mm Mains connection: 230 V EU plug

## Additional Information

---

No.	BL-0158
Manufacturer (Brand)	Bandelin
VGKL number	100733242
Transport temperature	Room temperature

