



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## **Bandelin® SONOREX DIGITEC DT 510 H-RC Ultrasonic bath with heating and infrared interface**

**€1,810.00  
plus VAT &  
Shipping**

### **Product Images**



## Description

---

High-performance ultrasonic cleaner with heating for use in the practice, laboratory and workshop. The proven and particularly robust SONOREX ultrasonic baths clean thoroughly and gently through intensive cavitation bubbles, even in the hard-to-reach places of the cleaning objects. A user-friendly and easy-to-wipe button operation is used for easy setting of the cleaning time and bath temperature. The remaining cleaning time as well as the current set and actual values of the bath temperature can be clearly read on the display. The fully automatic SweepTec function prevents standing waves in the ultrasonic field and achieves particularly uniform cleaning. The degassing function (Degas), which can be activated at the beginning of the cleaning process, additionally enhances the cleaning effect of the ultrasound. Suitable baskets, lids and cleaning and disinfection concentrates are available as accessories for optimally coordinated ultrasonic cleaning. Internal dimensions L x W x D: 300 x 240 x 150 mm Mains connection: 230 V EU plug

## Additional Information

---

No.	BL-0114
Manufacturer (Brand)	Bandelin
VGKL number	100733091
Transport temperature	Room temperature

