



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bandelin® SONOREX DIGITEC DT 514 BH-RC Ultrasonic bath with heating and infrared interface

€2,490.00
plus VAT &
Shipping

Product Images



Description

High-performance ultrasonic cleaner with heating for use in the practice, laboratory and workshop. The proven and particularly robust SONOREX ultrasonic baths clean thoroughly and gently through intensive cavitation bubbles, even in the hard-to-reach places of the cleaning objects. A user-friendly and easy-to-wipe button operation is used for easy setting of the cleaning time and bath temperature. The remaining cleaning time as well as the current set and actual values of the bath temperature can be clearly read on the display. The fully automatic SweepTec function prevents standing waves in the ultrasonic field and achieves particularly uniform cleaning. The degassing function (Degas), which can be activated at the beginning of the cleaning process, additionally enhances the cleaning effect of the ultrasound. Suitable baskets, lids and cleaning and disinfection concentrates are available as accessories for optimally coordinated ultrasonic cleaning. Internal dimensions L x W x D: 325 x 300 x 200 mm Mains connection: 230 V EU plug

Additional Information

No.	BL-0116
Manufacturer (Brand)	Bandelin
VGKL number	100733095
Transport temperature	Room temperature

