



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bandelin® SONOREX PR 140 DH Ultrasonic bath with heater

**€4,227.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The SONOREX pipette cleaner with heating gently and thoroughly removes dirt and incrustations by means of intensive ultrasonic cleaning. The device is also suitable for pipettes made of plastic, burettes and other glass devices. Features:

- Gentle pipette cleaning of glass volumetric instruments as well as parts up to a length of 755 mm
- Optional heating (DH) - thoroughly dissolves grease residues
- Particularly robust due to stainless steel tub made of 1.4301
- Tilting handles for easier transport
- Can be placed on the floor near the drain
- Cleaning liquid can be used several times, fast emptying due to ball valve.

Ready-to-use set: - SONOREX pipette cleaner PR 140 D - pipette basket K 140 B - lid D 140 D - cleaning concentrates TICKOPUR R 33 - 5 liters TICKOPUR TR 3 - 1 liter

Inner dimensions L x W x D: 150 x 150 x 895 mm  
Mains connection: 230 V EU plug

## Additional Information

---

No.	BL-0033
Manufacturer (Brand)	Bandelin
VGKL number	100732070
Transport temperature	Room temperature

