



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bandelin® SONOREX SUPER RK 1040 Ultrasonic bath 39.5 liters

**€3,719.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

High-performance ultrasonic cleaners for use in the practice, laboratory and workshop. The proven and particularly robust SONOREX ultrasonic baths clean thoroughly and gently through intensive cavitation bubbles, even in the hard-to-reach places of the cleaning objects. A user-friendly rotary handle control is used for easy adjustment of the cleaning time. The fully automatic SweepTec function prevents standing waves in the ultrasonic field and achieves particularly even cleaning. Matching baskets, lids and cleaning and disinfection concentrates are available as accessories for optimally coordinated ultrasonic cleaning. Internal dimensions L x W x D: 500 x 500 x 195 mm Mains connection: 230 V EU plug

## Additional Information

---

No.	BL-0125
Manufacturer (Brand)	Bandelin
VGKL number	960723014
Transport temperature	Room temperature

