



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bandelin® SONOREX SUPER RK 156 Ultrasonic bath 6 liters

€1,351.00
plus VAT &
Shipping

Product Images



Description

High-performance ultrasonic cleaners for use in the practice, laboratory and workshop. The proven and particularly robust SONOREX ultrasonic baths clean thoroughly and gently through intensive cavitation bubbles, even in the hard-to-reach places of the cleaning objects. A user-friendly rotary handle control is used for easy adjustment of the cleaning duration. Suitable baskets, lids and cleaning and disinfection concentrates are available as accessories for optimally coordinated ultrasonic cleaning. Internal dimensions L x W x D: 500 x 140 x 100 mm Mains connection: 230 V EU plug

Additional Information

No.	BL-0090
Manufacturer (Brand)	Bandelin
VGKL number	960723027
Transport temperature	Room temperature

