



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bandelin® SONOREX SUPER RK 514 BH ultrasonic bath with heater 18.7 liters

**€2,351.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

High-performance ultrasonic cleaner with heating for use in the practice, laboratory and workshop. The proven and particularly robust SONOREX ultrasonic baths clean thoroughly and gently through intensive cavitation bubbles, even in the hard-to-reach places of the cleaning objects. A user-friendly rotary handle control is used for easy adjustment of the cleaning duration and bath temperature. Matching baskets, lids and cleaning and disinfection concentrates are available as accessories for optimally coordinated ultrasonic cleaning. Internal dimensions L x W x D: 325 x 300 x 200 mm Mains connection: 230 V EU plug

## Additional Information

---

No.	BL-0045
Manufacturer (Brand)	Bandelin
VGKL number	960723010
Transport temperature	Room temperature

